

# GC%

- ✓ GC content for most eukaryotes: 40% - 50%
- ✓ GC content for bacterial genome: <15% - >75%

# Per Base Sequence Quality

- X-axis: position indexes
- Y-axis: Phred quality scores
- Blue line represents the mean quality scores
- The quality scores on the:
  - ✓ green region ( $Q > 28$ ) are very good
  - ✓ orange region ( $20 < Q < 28$ ) are reason-able
  - ✓ red region ( $Q < 20$ ) are poor

# Per Base Sequence Content

- X-axis: position indexes
- Y-axis: base percentage
- Higher percentage of some bases at the beginning of the x-axis may indicate contaminating remnants of adaptor sequences or other contaminating sequences.
- If the difference between any of the four bases in any position:
  - ✓ greater than 10%: a warning message
  - ✓ greater than 20%: failure of this metric

# Per Sequence GC Content

- X-axis: mean GC percentage per read
- Y-axis: number of reads
- A warning sign is displayed if the observed distribution deviates from normal distribution by a sum of more than 15% of the reads.
- A failure sign will be displayed if the distribution deviates by a sum of more than 30% of reads.

# Per Base N Content

- X-axis: position indexes
- Y-axis: N base percentages
- A warning is issued if any position shows an N content of greater than 5% and a failure sign if any position shows an N content of greater than 20%.

# Sequence Length Distribution

- X-axis: sequence length
- Y-axis: frequency
- A warning is displayed if the reads do not have the same length.

# Overrepresented Sequences

- A warning will be issued if a sequence is overrepresented more than 0.1% of the total
- A failure will occur if the overrepresentation is more than 1% of the total.

# Adapter Content

- X-axis: position indexes
- Y-axis: frequency
- A warning is raised if any sequence is present in more than 5% of all reads
- A failure occurs if any sequence is present in more than 10% of all reads.



# FastQC Command

- `fastqc filename.fastq.gz -o output_directory`

e.g., `fastqc SRR19551358.fastq.gz -o output/`

# Trimmomatic command

- TrimmomaticSE -threads 4  
SRR19551358.fastq.gz  
output/trimming/SRR19551358\_trimmed.fastq.g  
z ILLUMINACLIP:adapters.fa:2:30:10  
LEADING:28 TRAILING:28 HEADCROP:15  
CROP:172 SLIDINGWINDOW:4:15 MINLEN:36

# Hisat2 command

- `hisat2-build SRR19551358.fasta  
output/mapping/genome_index`
- `hisat2 -p 4 -x output/mapping/genome_index -U  
output/trimming/SRR19551358_trimmed.fastq.gz  
-S output/mapping/SRR19551358_aligned.sam --  
summary-file  
output/mapping/alignment_summary.txt`